

Sclavo-SED 10

Erythrocyte sedimentation rate analyser by a modified Westergren method



- ESR performed on the same Cell Blood Count (CBC) sample
- No need of dedicated Citrate tubes
- Innovative calculation algorithm supported by neural network
- Continuous loading of samples
- Results in only 25 minutes
- 10 inches Touch Screen user interface
- No production of liquid waste
- Safe and environmentally friendly



MADE IN ITALY

sclavo
diagnostics

Via Po 26/28, Loc. Pian dei Mori – 53018 Sovicille (Siena), Italy

Tel. +39 0577 39 041 | Email: sclavoexport@sclavo.it | servizio.clienti@sclavo.it | ordini@sclavo.it